



**MISSOURI DEPARTMENT OF TRANSPORTATION  
MATERIALS ENGINEERING  
Jefferson City, Missouri**

**Test Method  
MoDOT T52  
DETERMINATION OF MAGNESIUM CHLORIDE  
IN MAGNESIUM CHLORIDE BRINE**

**1.0 SCOPE**

**1.1** This method describes a procedure for determining the percent Magnesium Chloride in Magnesium Chloride brines by atomic absorption spectroscopy.

**2.0 REAGENTS AND APPARATUS**

**2.1** An Atomic Absorption Spectrophotometer.

**2.2** Hydrochloric Acid (HCl), Sp. Gr. 1.19.

**2.3** Lanthanum Oxide ( $\text{La}_2\text{O}_3$ ), reagent grade.

**2.4** 1000 ppm Magnesium stock solution. This solution can be purchased from a number of sources or it can be prepared in the laboratory from a suitably pure Magnesium salt.

**3.0 PREPARATION OF STOCK SOLUTIONS**

**3.1** Lanthanum stock solution - 10 percent. Add 200 ml of distilled water to 117.28 g of Lanthanum Oxide ( $\text{La}_2\text{O}_3$ ) in a 2000 ml beaker. While stirring, carefully add 500 ml of Hydrochloric Acid. When solution is complete, cool to room temperature, transfer to a 1000 ml volumetric flask and dilute to volume with distilled water.

**3.2** 30 ppm Magnesium stock solution. Pipette 30 ml of 1000 ppm Magnesium stock solution into a 1000 ml volumetric flask and dilute to volume with distilled water.

**4.0 PREPARATION OF STANDARD SOLUTIONS**

**4.1** Blank solution. Pipette 10 ml of Lanthanum stock solution into a 100 ml volumetric flask and dilute to volume with distilled water.

**4.2** 23.5%  $\text{MgCl}_2$ . Pipette 20 ml of 30 ppm Magnesium stock solution and 10 ml of Lanthanum stock solution into a 100 ml volumetric flask and dilute to volume with distilled water.



**4.3** 29.3%  $\text{MgCl}_2$ . Pipette 25 ml of 30 ppm Magnesium stock solution and 10 ml of Lanthanum stock solution into a 100 ml volumetric flask and dilute to volume with distilled water.

## **5.0 PROCEDURE**

**5.1** Weigh, to the nearest 0.1 mg, 0.5 g of brine into a 1000 ml volumetric flask and dilute to volume with distilled water. Transfer by pipette, a 20 ml aliquot to a 100 ml volumetric flask. Pipette 10 ml of Lanthanum stock solution into the 100 ml volumetric flask and dilute to volume with distilled water.

**5.2** Calibrate the instrument using the blank solution and the 23.5% and 29.3% Magnesium standard solutions, then determine %  $\text{MgCl}_2$  on the sample solution.

## **6.0 CALCULATION AND REPORT**

**6.1** The method of calculation will vary with the make and model of the instrument used. Report the percent Magnesium Chloride to the nearest 0.1% as follows:

% Magnesium Chloride ( $\text{MgCl}_2$ )

