



**MISSOURI DEPARTMENT OF TRANSPORTATION
MATERIALS ENGINEERING
Jefferson City, Missouri**

**Test Method
MoDOT T25
DETERMINATION OF CHLORINE IN
CHLORINATED RUBBER AND CHLORINATED PARAFFIN**

1.0 SCOPE

This method describes a procedure for determining the percent Chlorine in Chlorinated Rubber and Chlorinated Paraffin used in traffic paint vehicle.

2.0 REAGENTS AND APPARATUS

- (a) Parr Peroxide Bomb Apparatus, 22 mL, flame ignition Macro Bomb, including ignition shield.
- (b) Sodium Peroxide (Na_2O_2), Reagent Grade, granular, low chloride.
- (c) Nitric Acid (HNO_3), Reagent Grade.
- (d) Silver Nitrate (AgNO_3), Reagent Grade.
- (e) Sodium Chloride (NaCl), Reagent Grade, dried at 105-110 C for 1 hour prior to use.
- (f) Millivolt Meter equipped with a Combination Chloride Electrode.

3.0 PREPARATION OF STANDARD SOLUTIONS

- (a) Standard Sodium Chloride Solution (0.0100N) - Weigh 0.5844 gm dried NaCl , Reagent Grade, dissolve in water and dilute to 1L.
- (b) Standard Silver Nitrate Solution (0.01N) - Weigh 1.699 gm AgNO_3 , Reagent Grade, dissolve in water and dilute to 1L. Standardize against 0.0100 N NaCl .

4.0 PROCEDURE



- (a) Weigh, to the nearest 0.1 mg., 0.40 gm sample for Chlorinated Paraffin or 0.30 gm sample for Chlorinated Rubber and transfer to a Parr Bomb cup which contains 15 gm Na_2O_2 . Thoroughly mix the contents of the fusion cup with a stirring rod, wipe the rod with a small piece of filter paper and place the paper in the cup. Place the bomb assembly in the protective shield and fire strongly with a Bunsen burner for 20 minutes. Cool to room temperature and disassemble. (NOTE: If difficulties arise in obtaining complete combustion of Chlorinated Paraffin samples, small amounts of powdered sucrose or powdered NaNO_3 may be mixed with the Na_2O_2 and sample the fusion cup prior to ignition. However, the total combustible matter should not exceed 0.5 gm.) Rinse the bomb cover into a 600-mL beaker with about 150 mL of water, place the cup on its side in the beaker, cover with a watch glass and warm gently until solution of the fused material is complete. Remove the cup and rinse well with water. Add HNO_3 in small increments until the solution is acid to methyl orange. Add 5 mL HNO_3 in excess and warm gently until solution is complete. Filter through a coarse paper into a 1000-mL volumetric flask, washing well with hot water. Allow to cool to room temperature and dilute to volume with water.
- (b) Transfer a 50-mL aliquot to a 250-mL beaker and determine percent chlorine by potentiometric titration.
- (c) Perform a blank determination following the procedure given in (a) and (b) above, omitting only the addition of sample.

5.0 CALCULATION AND REPORT

Report percent chlorine to the nearest 0.1% as follows:

$$\% \text{ Chlorine (Cl)} = \frac{(A - B) \times C \times 0.709}{D} \times 100$$

Where:

- A = milliliters of AgNO_3 solution to titrate sample,
B = milliliters of AgNO_3 solution to titrate blank,
C = normality of AgNO_3 solution,
D = grams of the sample.

